




CAVITATION
Non-invasive body contouring

Cavitation is a selective non-invasive technology that uses low frequency ultrasound (30-70 KHz). It is used in aesthetic medicine for non-surgical reduction of localized fat and cellulite (P.E.F.S.). The cavitation, painless and free of side effects, can reduce outpatient volumes and circumferences with anesthesia and recovery time. The cavitation phenomenon develops in depth and, thanks to its lipolytic action associated with increased venous-lymphatic circulation, allows to obtain an improvement of the silhouette and the skin elasticity without damaging the surrounding tissues.



PURE:SKIN
 INTERNATIONAL

TECHNICAL CHARACTERISTICS

Display	6" colour touch screen
Adjustable treatment time	1-30 min
Emission Frequency	38 KHz
Stored protocols	11
Storable protocols	200
Output channels	2 independent - Ultracav 2100 1 - Ultracav Touch Mobile
Supplied probe	1 applicator 25cm2
Dimensions (LxWxH)	39 x 31 x 86 cm - Ultracav 2100 39 x 30 x 17cm - Ultracav Touch Mobile
Weight	28Kg - Ultracav 2100 3.5Kg - Ultracav Touch Mobile

PROVEN EFFICIENCY

Focussed applications in total safety

The cavitation equipment signed PureSkin International are the result of careful experimental work made by our scientific consultants and represent the achievement of an important goal: to offer a non-invasive technology against cellulite (P.E.F.S.) and localised fat.



Before



After



Before



After



PURE:SKIN
 INTERNATIONAL


Ultra Cav

Medical grade cavitation

Ultrasound against localised fat and cellulite (P.E.F.S.)



MADE IN ITALY



TARGETED ACTION FOR IMMEDIATE EFFECTS

- Lipoclastic action and progressive remodeling of the silhouette
- Reduction of orange peel appearance and elimination of fibro-fatty nodules of cellulitis
- Body contouring
- Oxygenation and revascularization of asphyxiated skin
- Drainage of stagnation liquids
- Reactivation of the peripheral circulation
- Absorption of revitalising or lipolytic actives
- Improvement of the tone and elasticity of tissues



GUARANTEED RESULTS

Ultracav 2100 and Ultracav Touch Mobile reduce cellulite blemishes (P.E.F.S.) and localised fat shaping the body in a painless and non-invasive way in contrast with other methods like aesthetic liposuction.



EI10101 - ULTRACAV TOUCH MOBILE

Portable equipment for cavitation
Certified CE 0476 for medical use

